

K2 An Analysis of Cannabinoids in Whole Blood

Adrian M. Taylor, PhD, SCIEX, 71 Four Valley Drive, Concord, ON L4K4V8, CANADA*

After attending this presentation, attendees will better understand how to effectively use protein precipitation and 2D Ultra High-Performance Liquid Chromatography/Mass Spectrometry (2D-UHPLC/MS) for the analysis of Δ^9 -Tetrahydrocannabinol (THC), 11-OH- Δ^9 -Tetrahydrocannabinol (THC-OH), and 11-nor-9-carboxy- Δ^9 -Tetrahydrocannabinol (THC-COOH) in whole blood.

This presentation will impact the forensic science community by illustrating how combining a very simple sample preparation with a 2D-LC analytic approach in UHPLC allows for reaching the required sensitivity cut-offs with high-throughput capacity and robustness. This method claims to be a very interesting substitute to a classic technique based on Gas Chromatography/Mass Spectrometry (GC/MS) for the quantitation on cannabinoids in whole blood.

Hypothesis: Protein precipitation and 2D-UHPLC-MS provides effective clean-up of whole blood to analyze THC, THC-OH, and THC-COOH in a semi-automated way.

Methods: Sample Preparation – 100 μ L of whole blood (or calibrators, quality controls and unknowns) are mixed with 20 μ L of internal standard mixture in methanol (containing deuterated analytes). The protein precipitation step is done by the addition of 200 μ L of cold acetonitrile, 0.1% formic acid. After a ten minute incubation in the freezer, samples were centrifuged for five minutes at 13,000rpm. Next 200 μ L of supernatant were mixed with 200 μ L of water (0.1% formic acid) and were analyzed. The injection volume for each sample was defined at 100 μ L. **UHPLC Conditions** – Chromatographic separation was achieved using two Kinetex™ columns (PFP, 5 μ m and XB-C18, 5 μ m). A Sciex 4500 QTRAP® was used for the MRM detection with an electrospray interface in positive/negative ionization mode. The run time of 2D-UHPLC program with equilibration was optimized at ten minutes using a solvent gradient of water 0.1% formic acid and methanol 0.1% formic acid.

Results: The sample preparation consists of a simple pre-treatment step using precipitation of whole blood samples, followed by dilution and injection. UHPLC methodology consists of a water/methanol mobile phase with an efficient 2D-LC gradient on Kinetex™ columns. Limits of quantitation were obtained below 0.5ng/mL for THC and its metabolites. Real samples were analyzed by this workflow. Results from a negative and a positive sample will be shown. Results are based on spiked standard curves in drug-free whole blood. The calculated concentrations for the positive sample were: 1.4ng/mL for THC, <0.5ng/ml for THC-OH, and 3ng/mL for THC-COOH.

Conclusion: Cannabinoids in whole blood analysis has been demonstrated on a 2D-UHPLC system, coupled to the SCIEX 4500 QTRAP® Mass spectrometer. Analysis is sensitive to detect concentration below 0.5ng/mL for each compound with Signal to Noise (S/N) greater than 10:1 (based on 3SD of the noise).

This method uses only 100 μ L of whole blood without offline solid phase extraction or liquid-liquid extraction. The short runtime allows for a high throughput whilst maintaining chromatographic performance batch after batch. No cross contamination was observed between the highest sample concentrations to the negative samples.

2D-UHPLC/MS, THC, Whole Blood