

## B52 A Single Assay for Human-Specific Quantification of Less Than One Picogram DNA and Detection of the Presence of PCR Inhibitors in Forensic Samples

Maura T. Costello, MFS, and James W. Schumm, PhD\*, The Bode Technology Group, 7364 Steel Mill Drive, Springfield, VA 22150

After attending this presentation, attendees will gain an improved method for quantification of human DNA that will detect less than one genome.

This presentation will impact the forensic community and/or humanity by demonstrating the how the application of the method described will lead to increased success in DNA analysis of samples con- taining very little DNA. Furthermore, the method will facilitate detection of sample components that inhibit amplification. This combination will provide critical information on first use of the sample (*e.g.*, dilution vs. using more sample), an especially critical choice when evaluating samples containing very little material.

This presentation describes the development, validation, and appli- cation of a duplex real-time PCR assay for human-specific quantification of DNA samples containing as little as 0.5 pg/µl of DNA. The assay simul- taneously detects PCR inhibitors within the sample.

It is important to include human-specific quantification of DNA in casework sample analysis to insure successful DNA amplification and pro- filing. Much recent research has focused on the use of real-time quantitative PCR to achieve this goal. This approach is less labor intensive, less time consuming, more accurate, and lends itself to automation better than previous methods such as slot-blot hybridization (1). This work builds on that described by Nicklas and Buel (2), Richard et al. (3), and the commer- cially available Quantifiler<sup>™</sup> Kit (Applied Biosystems, Foster City, CA). Researchers combined the sensitivity and human specificity of *Alu*-based real-time quantification with the presence of an internal positive control allowing detection of PCR inhibitors in the sample.

*Alu* sequences are short, repeated elements that are interspersed throughout the primate genome in upwards of 500,000 copies. The Yb8 subfamily of *Alu* genes was selected because of its sequence specificity to higher primates (4). Using this target, primers and a fluorogenic probe were developed for a quantitative real-time PCR assay (5). The assay also contains an internal positive control (IPC) system that is multiplexed with the *Alu* quantification system, consisting of a fixed quantity of non-human DNA template added to each reaction well, and a second set of primers and fluorogenic probe specific for the non-human template. The combination of human DNA quantity data from the *Alu* system and DNA quality data from the IPC system provides the analyst with substantial information to aid in deciding dilution or concentration schemes prior to STR amplification, thereby significantly reducing the number of samples that need to be re-evaluated following initial profiling.

Validation work indicates the assay is accurate and precise in the range of 50 ng/µl to 0.5 pg/µl. Thus less than one human genome equivalent can be detected accurately. Species specificity tests indicate the assay is at least 5000 times more specific for higher primate DNA than any other species tested. The IPC system is very sensitive to inhibition observed with addition of hematin, indigo, or humic acid. The assay has been successful with a variety of non-probative sample types.

The features of this assay will allow us to apply it very effectively to evaluation of touch evidence samples. With so little sample available in these situations, it is critical to make the right decision to use more (with limiting amounts) or less (with inhibitors present) extracted DNA in the first profiling test.

## References:

- 1. Walsh, P.S.; Varlaro, J.; Reynolds, R.A.; Rapid Chemiluminescent Method for Quantitation of Human DNA. *Nucleic Acids Res.* 1992; 20:5061-5.
- Nicklas, J.A.; Buel, E.; Development of an ALU-based, Real-time PCR Method for Quantitation of Human DNA in Forensic Samples. *J Forensic Sci.* 2003; 48(5):935-44.
- 3. Richard, M.L.; Frappier, R.H.; Newman, J.C.; Developmental Validation of Real-Time Quantitative PCR Assay for Automated Quantification of Human DNA. *J Forensic Sci.* 2003; 48(5):1041-6.
- 4. Carroll, M.L.; et al. Large-Scale Analysis of the Alu Ya5 and Yb8 Subfamilies and Their Contribution to Human Genomic Diversity. *J Mol Biol.* 2001; 311:17-40.
- Holland PM, Abramson RD, Watson R, Gelfand, DH. Detection of Specific Polymerase Chain Reaction Product by Utilizing the 5'-3' Exonuclease Activity of Thermus Aquaticus DNA Polymerase. Proc Natl Acad Sci USA. 1991; 88:7276–7280.

Quantification, Low Copy Number, Inhibition