



A108 Raman Microspectroscopy and Advanced Statistics for Detection and Characterization of Gunshot Residue (GSR)

Igor K. Lednev, PhD*, Justin Bueno, MS, and Vitali Sikirzhyski, MS, SUNY Albany, 1400 Washington Ave, Albany, NY 12222

The goal of this presentation is to describe the development of a novel and alternative method for Gunshot Residue (GSR) analysis.¹ After attending this presentation, attendees will have a better understanding of recent advancements of this application of Raman microspectroscopy for GSR analysis, identification, and discrimination. The implementation of advanced statistics to differentiate experimental Raman spectra collected from non-equivalent GSR samples will be discussed.

This presentation will impact the forensic science community by having the potential to greatly impact the accuracy and effectiveness of shooting incident investigations.

Raman spectroscopy has numerous applications in forensic chemistry. Raman analysis is a technique which can obtain confirmatory class identification of analytes through low intensity laser light scattering. The technique is non-destructive, rapid, sensitive, and requires little or no sample preparation. Furthermore, portable Raman spectrometers are readily available, allowing for crime scene accessibility. Raman spectroscopy offers several advantages over the current methodology for GSR analysis. The technique has been shown to detect components from both the organic and inorganic constituents of GSR. This is contrary to current GSR elemental analysis methods which rely solely on the detection of the heavy metals (lead, barium, and antimony). This is problematic since environmental concerns have led to the increased popularity in heavy metal-free or "green" ammunition. It has been found that in the absence of heavy metals, current elemental analysis techniques are severely hindered when making accurate identification of GSR samples. Additionally, the probability of environmental and manufacturing particles assigned (incorrectly) as being GSR has increased with the onset of "green" ammunition. Until recently, the application of Raman spectroscopy for GSR analysis was largely unexplored, although this approach is not dependent upon detecting metals, and is more capable of differentiating environmental contaminants and GSR. Therefore, a Raman spectroscopic method displays numerous advantages in specificity when compared to current techniques.

The firearm discharge process could be considered analogous to a complex chemical reaction. Therefore, the chemical composition of the products (GSR particles) is directly related to the chemical nature of the reagents (firearm-ammunition combination) and the conditions of the reaction. Preliminary results show that Raman data collected from GSR particles originating from different firearm-ammunition discharges were successfully classified according to caliber; 0.38 inch and 9mm caliber firearm discharge samples were probed using a 785-nm Raman excitation. Resulting data was treated with statistical methods (performed using Matlab with the PLS toolbox) such as Principle Component Analysis (PCA) and Support Vector Machines (SVM). The results show a high probability of this method to correctly classify data from the two examined calibers. Preliminary results illustrate that the variations between non-equivalent GSR samples can be detected through this method. Since GSR is often collected from a suspect, the application of this method to forensic investigations would provide a link between GSR collected from the shooter and the crime scene.

This emerging technique illustrates the possibility for an on-scene, non-destructive, identification and chemical characterization method for GSR. This method has the potential to greatly impact the forensic science community by increasing the accuracy (and discriminatory power) of GSR detection. The most direct application for this research is a method to exclude a specific firearm-ammunition combination as producing an evidentiary GSR sample. The comparison of a laboratory-generated GSR sample discharge and an evidentiary GSR sample, can be made without extensive preliminary studies.

Reference:

- ¹ Bueno J, Sikirzhyski V, Lednev IK. Raman spectroscopic analysis of gunshot residue offering great potential for caliber differentiation. *Anal Chem* 2012;84:4334-4339.

Gunshot Residue, Spectroscopy, Statistical Analysis